

Protocol submission information for EQA of NordiQC

When filling out the protocol forms for the NordiQC EQA's, please use below information to fill the different lines.

Please check on your Genie instrument if the antibody in question is using Citrate or high pH buffer as pretreatment.

New protocol - 73 - ALK (lung)

Staining platform	
Staining platform	Sakura Tissue-Tek Genie ▼

Primary antibody	
Primary antibody clone	Sakura Finetek (8330-C010) - 1A4 ▼
Lot number	8330..... fill in a lot number you have on stock
Incubation time (minutes)	30
Incubation temperature (Celcius)	20

Epitope Retrieval, HIER	
Epitope retrieval, HIER	<input checked="" type="radio"/> YES <input type="radio"/> NO
Device	Other ▼
HIER buffer	Sakura - Tissue-Tek Genie High pH Antigen Retriev ▼
Efficient Heating Time (minutes)	45
Max. heating temperature (Celcius)	98

Epitope Retrieval, proteolysis	
Epitope retrieval, proteolysis	<input type="radio"/> YES <input checked="" type="radio"/> NO

Visualization system	
Visualization system	Tissue-Tek Genie Pro Detection Kit, DAB - 8826-K ▼
Incubation time linker (minutes)	15
Incubation time polymer (minutes)	20
Incubation temperature (Celcius)	20

Chromogene	
Chromogene	Sakura - Tissue-Tek Genie Pro Detection Kit, DAB I ▼
Incubation time (minutes)	6
Incubation temperature (Celcius)	25
Enhancement	CuSO4 ▼

Comment	
Comments to protocol	<input type="text"/>