

# Tissue-Tek Genie®

## Dewax Solution

### Instructions for use

#### Intended use

For *in vitro* diagnostic use.

Tissue-Tek Genie® Dewax Solution is designed to remove paraffin wax from paraffin embedded specimen sections on the Tissue-Tek Genie® Advanced Staining System. The clinical interpretation must be made in conjunction with histological examination, relevant clinical information, other diagnostic tests and proper controls by a qualified pathologist.

#### Limitations

This product has been optimized on the Tissue-Tek Genie Advanced Staining System, using Tissue-Tek Genie® reagents and FFPE specimen sections. Staining quality may diminish when used with other systems and/or reagents.

#### Summary and principle

FFPE specimen sections mounted on microscope slides are treated with Tissue-Tek Genie Dewax Solution on the Tissue-Tek Genie Advanced Staining System prior to the immunohistochemistry and/or *in situ* hybridization staining processes. Specimen sections in the presence of Tissue-Tek Genie Dewax Solution are gently heated for optimized removal of paraffin.

#### Cautions and warnings

For professional use only. Take reasonable precautions when handling. Avoid contact of reagents with eyes,

skin, and mucous membranes. Wear protective gloves, clothing, and eye/face protection.

Bottles filled with reagent are for multiple use. Do not attempt to re-fill or add additional reagent. Discard the bottle when empty.

It is recommended to include appropriate controls on each specimen slide to help in identifying any deviation that might occur during the staining process.

All disposal practices must be in compliance with all Federal, State/Provincial and local laws and regulations. Refer to the SDS for further information.

#### Storage conditions

Store this product at 15-30°C.

#### Instructions for use

1. Remove one bottle and one draw tube from the box and the bag.
2. Scan the 2D barcode on the bottle label to register the lot with the software.
3. Remove the black cap from the bottle.
4. Place the draw tube into the bottle and tighten the cap.
5. Place the Tissue-Tek Genie Dewax Solution line (line one of the Tissue-Tek Genie Advanced Staining System) into the bottle and tighten the black connector attached to the line. Ensure the line reaches the bottom of the bottle.

6. Gently place the bottle onto the Tissue-Tek Genie Dewax Solution slot of the Tissue-Tek Genie Advanced Staining System.
7. Prime the Tissue-Tek Genie Dewax Solution line if necessary.

### Material required but not supplied

The following reagents may be required for staining but are not provided:

- Tissue-Tek Genie® Citrate Antigen Retrieval Solution (REF 8742-G001)
- Tissue-Tek Genie® High pH Antigen Retrieval Solution (REF 8744-G001)
- Tissue-Tek Genie® Wash Buffer Solution (REF 8874-G004)
- Tissue-Tek Genie® Hematoxylin (REF 8830-M250)
- Tissue-Tek Genie® Non-Immune Mouse Ig Antibody, Negative Control (REF 8604-C010, 8604-M250)
- Tissue-Tek Genie® Non-Immune Rabbit Ig Antibody, Negative Control (REF 8605-C010, 8605-M250)
- Tissue-Tek Genie® DUO Non-immune Mouse and Rabbit Ig Antibody Cocktail, Negative Control (REF 8482-C010, 8482-M250)
- Tissue-Tek Genie® detection systems
- Tissue-Tek Genie® primary antibodies

Further information can be found on the Sakura Finetek USA website at [www.sakuraus.com/Genie](http://www.sakuraus.com/Genie)

### Order information

#### Product code, product name and quantity

REF 8865-G001 Tissue-Tek Genie® Dewax Solution, Ready-To-Use, 3.8 liters, 1 bottle; 1 unit.

**NOTE:** The Safety Data Sheet (SDS) is available online on the Sakura Finetek USA website at [www.sakuraus.com/SDS.html](http://www.sakuraus.com/SDS.html)

### Contact


If located within the United States, contact Sakura Finetek USA, Inc. by calling toll free **1-800-725-8723** or contact your Sakura Finetek representative or authorized distributor.

In countries, other than the United States, contact the nearest authorized Sakura Finetek instrument distributor or representative. Contact details may be found at [www.sakura.com](http://www.sakura.com)



## Symbols

|   |   |
|---|---|
|  | Catalog number                                      |
|  | Batch code  |
|  | <i>in vitro</i> diagnostic medical device           |
|  | Temperature limitation                              |
|  | Use by  |
|  | Manufacturer  |
|  | Consult instructions for use                        |
|  | European Conformity                                 |
|  | Authorized representative in the European Community |

Storage: 15°C  30°C







|   |  |
|---|--|
|  | Sakura Finetek USA, Inc.<br>1750 W 214 <sup>th</sup> Street<br>Torrance, CA 90501 U.S.A.       |
|   | Sakura Finetek Europe B.V.<br>Flemingweg 10a<br>2408 AV Alphen aan den Rijn<br>The Netherlands |
| Made in U.S.A.  |  |

GS-31840 Rev. A