

Tissue-Tek Genie®

anti-CD10 Rabbit Monoclonal Antibody [QR021]

REF 8386-C010

Instructions for use

For *in vitro* diagnostic use.

Intended purpose

Intended use: The Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021] is designed to qualitatively detect CD10 protein in formalin-fixed, paraffin-embedded (FFPE) human tissue specimen sections by immunohistochemistry (IHC) staining on the automated Tissue-Tek Genie® Advanced Staining System.

This device functions as an aid for diagnosis and shall be used by a qualified pathologist with a panel of other antibodies to classify a subset of B-cell lymphomas such as diffuse large B-cell lymphoma (DLBCL), follicular lymphoma and Burkitt lymphoma, and renal cell carcinoma.

Limitations

The Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021] has been optimized for use with the Tissue-Tek Genie® Advanced Staining System, Tissue-Tek Genie® reagents, and FFPE specimen sections. Staining quality may be diminished when used with other systems and/or reagents.

The clinical interpretation must be made in conjunction with histological examination, relevant clinical information, a panel of other antibodies, other diagnostic tests, and proper controls by a qualified pathologist.

Staining quality may be diminished by improper or incomplete removal of the paraffin.

The sensitivity of this antibody to identify the CD10 protein may be affected by improper specimen handling. This may alter antigenicity, weaken detection, and may generate false negative results.

Special processing of tissues such as decalcification of bone marrow tissues may lead to inconsistent staining.

Positively charged slides are recommended to obtain optimal staining with the Tissue-Tek Genie Advanced Staining System.

Summary and principle

Immunohistochemistry (IHC) staining is an established *in vitro* diagnostic method to visualize the presence of specific proteins expressed within a tissue section to study the microscopic features.

IHC staining is accomplished in two steps:

1) a primary antibody recognizes a particular target protein expressed on a specific cell compartment of a specific cell type on various tissues, and

2) a secondary and tertiary antibody conjugated to a chromogenic enzyme bind with the primary antibody for the detection of the antibody-antigen interaction. In chromogenic detection under a light microscope, an enzyme conjugated to the antibody cleaves a substrate to produce a colored precipitate at the location of the protein.

CD10, also known as common acute lymphoblastic leukemia antigen (CALLA), is a cell surface neutral endopeptidase. It exhibits a predominantly membranous and occasionally cytoplasmic staining pattern in a wide variety of hematopoietic and nonhematopoietic tissues.¹⁻⁷

In normal tissues, CD10 is expressed on the cell surface of bone marrow stem cells, a small subset of bone marrow immature B-cells and lymph node germinal centers, a subpopulation of parafollicular T-cells, and mature neutrophils.^{2,5,8-10} CD10 is also expressed in various non-lymphoid cells such as epithelial cells of glomeruli and proximal tubules of kidney, bile canaliculi of liver, alveolar epithelial cells of lung, ductal cells of prostate, trophoblastic cells of placenta, myoepithelial cells of secretory glands, endometrial stromal cells of uterus, and at the apical surface (brush border) of enterocytes in the small intestine.^{2-4,6,8,11}

In tumor tissues, CD10 expression has been detected in various types of neoplasms including in 44-100% of germinal center B-cell (GCB) diffuse large B-cell lymphoma (DLBCL), 0-5% of non-GCB DLBCL,¹²⁻¹⁴ 46-100% of follicular lymphoma (FL),^{15,16} 79-100% of Burkitt lymphoma (BL),¹⁷ 88-100% of renal cell carcinoma (RCC) clear cell type,^{4,7,18} endometrial stromal sarcoma,⁴ hepatocellular

carcinoma (HCC),^{19,20} urothelial carcinoma,³ and prostate carcinoma.²¹

CD10 is a marker of B-cell differentiation and germinal center origin. The expression level of CD10 is relatively high in precursor B lymphoma/leukemias and some non-hematopoietic tumors.^{1,5,22} Anti-CD10 immunohistochemistry is useful in characterizing a subset of malignant lymphomas derived from germinal center B-cells (GCB) such as DLBCL, FL, and BL. It is also useful in identifying RCC, including the clear cell and papillary variants.^{1,5,7,22}

Anti-CD10 antibody labels CD10 in both normal and neoplastic cells and has a predominantly membranous and occasionally cytoplasmic staining pattern.

The Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021] is a primary antibody against the human CD10 protein and is provided in buffered saline containing 1% bovine serum albumin and 0.09% sodium azide. FFPE specimen sections are placed on positively charged slides and the paraffin is removed using the Tissue-Tek Genie® Dewax Solution (REF 8865-G001), after which heat-induced epitope retrieval is performed using the Tissue-Tek Genie® High pH Antigen Retrieval Solution (REF 8744-G001).

IHC demonstration of CD10 protein in FFPE specimen sections is achieved through the use of the Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021] and the Tissue-Tek Genie® PRO Detection Kit, DAB (REF 8826-K250). This procedure entails the sequential application of antibody and kit components as follows:

- Tissue-Tek Genie® Protein Block
- Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021]
- Tissue-Tek Genie® Peroxidase Block
- Tissue-Tek Genie® Link (binds to the primary antibody)
- Tissue-Tek Genie® Poly HRP-Conjugate (binds to the Link)
- Tissue-Tek Genie® DAB (visualizes the detected protein)
- Tissue-Tek Genie® DAB Intensifier

Tissue-Tek Genie® Hematoxylin (REF 8830-M250) is then used to visualize the nuclei of cells. The IHC stained slide is coverslipped and reviewed using a light microscope.

Expected results

Cellular staining pattern: membranous and cytoplasmic

Positive control tissue: tonsil, appendix, liver, kidney, selective B-cell lymphoma

The specificity and intended use of this antibody were validated by performing IHC staining on the Tissue-Tek Genie Advanced Staining System using normal and neoplastic FFPE tissue sections as follows.

Analytical sensitivity/specificity: A total of 32 types and 128 specimens of normal FFPE tissues were tested. Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021] exhibited moderate to strong, predominantly membranous and cytoplasmic staining in various tissues such as germinal center B-cells in tonsil (17/17) and

appendix (11/11), bile canaliculi in liver (9/9), and proximal tubules and glomeruli in kidney (14/14). No staining was observed in mantle zone B-cells and squamous epithelial cells in tonsil (0/17). Positive staining was observed in alveolar epithelial cells of lung (4/4), glandular cells of prostate (5/5), trophoblastic cells of placenta (5/5), the apical surface (brush border) of enterocytes of the small intestine (2/2), stromal cells of uterus (1/2), myoepithelial cells of secretory glands (7/7) including mammary, prostate, and salivary glands. Positive staining in endothelial cells and lymphocytes has also been observed.

Precision studies for Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021] lots were completed. FFPE tissue sections of tonsils were used. Studies were conducted to demonstrate reproducibility in lot-to-lot (minimum of 2 lots), run-to-run (minimum of 2 Genie runs), instrument-to-instrument (2 Genies), station-to-station (minimum of 2 stations), and operator-to-operator (2 operators). The results were compared and met their acceptance criteria: moderate to strong membranous staining in germinal center B-cells in tonsil.

These results demonstrate precision of the Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021], which was consistent across lots, runs, instruments, stations, and operators.

Diagnostic sensitivity/specificity: A total of 40 types and 150 specimens of neoplastic FFPE tissues were tested. Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021] demonstrated membranous positivity of neoplastic cells in DLBCL of GCB subtype (13/14, 92.8%), FL (8/8, 100%), BL (3/3, 100%), and RCC clear cell types (11/11, 100%). No staining was observed in carcinomas of the breast (0/10), ovary (0/3), thyroid (0/4), stomach (0/2), and Hodgkin lymphoma (0/3).

Diagnostic specificity was demonstrated by the absence of staining in neoplastic cells of DLBCL of non-GCB subtype (18/18, 100%).

Clinical performance:

The Tissue-Tek Genie anti-CD10 Rabbit Monoclonal Antibody [QR021] demonstrated conformity with the expected clinical performance through analytical studies and assessments of diagnostic performance in conjunction with established scientific validity (summarized in the “Summary and principle” section above with data from the references at the end of this IFU) based on information on other IVD medical devices with the same antibody, textbooks, and available peer-reviewed literature.

Tissue specimen	Established validity	Tissue specimen tested
Tonsil	CD10 is expressed in virtually all germinal center B-cells. ⁹	Moderate to strong, predominantly membranous, staining in germinal center B-cells (17/17) was observed.
Liver	CD10 staining is observed in bile canaliculi. ^{4,8}	Predominantly strong membranous and cytoplasmic staining of



		bile canaliculi was observed in liver (9/9).
Kidney	CD10 is expressed in proximal tubular and glomerular epithelial cells. ^{4,7,8}	Moderate to strong membranous staining was observed in proximal tubules and glomeruli in kidney (14/14).
DLBCL (GCB subtype)	CD10 is detected in 44-100% of neoplastic B-cells in DLBCL (GCB subtype) cases. ¹²⁻¹⁴	Moderate to strong membranous staining was observed in the neoplastic cells of DLBCL (GCB subtype) cases (13/14, 93%).
DLBCL (non-GCB subtype)	CD10 is negative in other neoplasms, such as DLBCL (non-GCB subtype) (0-5%) cases. ¹²⁻¹⁴	Absence of staining was observed in the neoplastic cells of DLBCL (non-GCB subtype) cases (18/18, 100%).
Follicular Lymphoma (FL)	CD10 is detected in 46-100% of neoplastic cells of follicular lymphoma cases. ^{15,16}	Moderate or strong membranous staining was observed in the neoplastic cells of follicular lymphoma cases (8/8, 100%).
Burkitt Lymphoma (BL)	CD10 is detected in 79-100% of neoplastic cells of Burkitt lymphoma cases. ¹⁷	Moderate or strong membranous staining was observed in the neoplastic cells of Burkitt lymphoma cases (3/3, 100%).
RCC (clear cell type)	CD10 is detected in 88-100% of neoplastic cells of renal cell carcinoma cases. ^{1,5,7,22}	Moderate or strong membranous staining was observed in the neoplastic cells of renal cell carcinoma cases (11/11, 100%).

Together, this information is sufficient to demonstrate conformity with relevant essential principles without the need for additional clinical performance data.

Cautions and warnings

For professional use only.

Take reasonable precautions when handling. Avoid contact of reagents with eyes, skin, and mucous membranes. Wear protective gloves, protective clothing, and eye protection.

The product should not be exposed to temperatures outside of the storage conditions.

Capsules filled with ready-to-use, pre-diluted, antibody are for single use only. The in-use stability of each capsule is approximately 72 hours when used outside the storage conditions. Do not attempt to refill or add additional reagent. Discard capsule after use.

It is recommended to include appropriate controls on each specimen slide to help in identifying any deviation that might occur during the staining process.

All disposal practices must be in compliance with all Federal, State/Provincial and local laws and regulations. Refer to the SDS for further information.

Specimen collection and preparation for analysis

Routinely processed, formalin-fixed, paraffin-embedded tissues are suitable for use with this reagent when used with Tissue-Tek Genie reagents and a Tissue-Tek Genie Advanced Staining System (see section "Material required but not supplied"). The recommended fixation is performed using 10% (v/v) neutral buffered formalin for 24-72 hours.²³ Variable results may occur because of prolonged fixation or special processes such as decalcification of bone marrow preparations. Each cut section should be 3-5 µm in thickness and placed on a positively charged glass slide. Slides containing the tissue section may be baked for at least 30 minutes to overnight (typically up to 16 hours) in a 58-60°C oven.²³

Storage conditions

Store this product at 2-8°C. Do not freeze. Return unused capsules to 2-8°C.

For the date of expiration, refer to the label on the product.

The reagent will be stable until its expiration date when stored and handled properly. Do not use the reagent beyond its assigned expiration date. Storage conditions other than those specified above must be verified by the user.

Do not use when precipitate is visible in the reagent container.

Instructions for use

Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021], RTU, 10 capsules/pack (REF 8386-C010):

1. Place the Tissue-Tek Genie® Reagent Dispensing Area Tag (RDA-Tag) attached to the capsule into the RDA.
2. Push the capsule into the RDA with foil side down and click the attached RDA-Tag down into place on the RDA.
3. Place the RDA on the desired station of the Tissue-Tek Genie Advanced Staining System.
4. Place the slide with the specimen section on the same station, specimen section side down.
5. Assign protocol 8386 to the same station.
6. Initiate execution of protocol 8386.
7. The RDA-Tag 8386 will be scanned and registered automatically when the staining process is initiated.
8. During the primary antibody application step, the antibody will be released from the capsule into the RDA and onto the specimen section on the slide.
9. The staining protocol continues to the end.

Material required but not supplied

- Tissue-Tek Genie® Advanced Staining System (REF 8200)
- Positive and negative tissue controls
- Drying oven capable of maintaining a temperature of 58-60°C
- Tissue-Tek Genie® Dewax Solution (REF 8865-G001)
- Tissue-Tek Genie® Wash Buffer Solution (REF 8874-G004)
- Tissue-Tek Genie® High pH Antigen Retrieval Solution (REF 8744-G001)
- Tissue-Tek Genie® Non-Immune Rabbit Ig Antibody, Negative Control (REF 8605-C010)
- Tissue-Tek Genie® PRO Detection Kit, DAB (REF 8826-K250)
- Tissue-Tek Genie® Hematoxylin (REF 8830-M250)
- Tissue-Tek Genie® Reagent Dispense Area [RDA] (REF 8616-G090)

Further information can be found on the Sakura Finetek USA website at www.sakuraus.com/Genie

Troubleshooting

Testing run should include proper reagent and tissue controls.

- If the positive control exhibits negative, or weaker, or stronger staining, or more background staining than expected, other positive controls on the same instrument run should be examined to determine if this is due to the antibody, other reagents, software, instrumentation, or the processing and fixation of tissue specimen(s).
- If the paraffin has not been removed completely, the deparaffinization procedure should be verified.
- If tissue sections have washed off, slides should be examined to ensure they are positively charged, and the specimen should be examined for possible inadequate processing or fixation.
- Refer to the Tissue-Tek Genie Advanced Staining System operating manual or contact your Sakura Finetek Technical support representative for information or assistance.

Order information / product provided

Product code, product name and quantity

REF 8386-C010 Tissue-Tek Genie® anti-CD10 Rabbit Monoclonal Antibody [QR021]; RTU, 10 capsules/pack.

NOTE: The Safety Data Sheet (SDS) is available online on the Sakura Finetek USA website at www.sakuraus.com/SDS.html

The Summary of Safety and Performance (SSP) is available online via EUDAMED.

References

1. Maguer-Satta V, Besancon R, Bachelard-Cascales E. Concise review: neutral endopeptidase (CD10): a multifaceted

environment actor in stem cells, physiological mechanisms, and cancer. *Stem Cells*. 2011; 29(3): 389-96. <https://doi.org/10.1002/stem.592>.

2. O'Malley DP, Fedoriw Y, Grimm KE, Bhargava P, Immunohistology of Lymph Node and Lymph Node Neoplasms, in *Diagnostic Immunohistochemistry*, D.J. Dabbs, Editor. 2021, Elsevier Health Sciences. p. 156-195.
3. Bahadir B, Behzatoglu K, Bektas S, Bozkurt ER, Ozdamar SO. CD10 expression in urothelial carcinoma of the bladder. *Diagn Pathol*. 2009; 4: 38. <https://doi.org/10.1186/1746-1596-4-38>.
4. Chu P, Arber DA. Paraffin-section detection of CD10 in 505 nonhematopoietic neoplasms. Frequent expression in renal cell carcinoma and endometrial stromal sarcoma. *Am J Clin Pathol*. 2000; 113(3): 374-82. <https://doi.org/10.1309/8VAV-J2FU-8CU9-EK18>.
5. Mishra D, Singh S, Narayan G. Role of B Cell Development Marker CD10 in Cancer Progression and Prognosis. *Mol Biol Int*. 2016; 2016: 4328697. <https://doi.org/10.1155/2016/4328697>.
6. Yasuda M, Itoh J, Satoh Y, Kumaki N, Tsukinoki K, Ogane N, Osamura RY. Availability of CD10 as a histopathological diagnostic marker. *Acta Histochemica et Cytochemica*. 2005; 38(1): 17-24.
7. Truong LD, Shen SS. Immunohistochemical diagnosis of renal neoplasms. *Arch Pathol Lab Med*. 2011; 135(1): 92-109. <https://doi.org/10.5858/2010-0478-RAR.1>.
8. McIntosh GG, Lodge AJ, Watson P, Hall AG, Wood K, Anderson JJ, Angus B, Horne CH, Milton ID. NCL-CD10-270: a new monoclonal antibody recognizing CD10 in paraffin-embedded tissue. *Am J Pathol*. 1999; 154(1): 77-82. [https://doi.org/10.1016/S0002-9440\(10\)65253-4](https://doi.org/10.1016/S0002-9440(10)65253-4).
9. Shipp MA, Stefano GB, Switzer SN, Griffin JD, Reinherz EL. CD10 (CALLA)/neutral endopeptidase 24.11 modulates inflammatory peptide-induced changes in neutrophil morphology, migration, and adhesion proteins and is itself regulated by neutrophil activation. *Blood*. 1991; 78(7): 1834-41.
10. Huang X, He C, Lin G, Lu L, Xing K, Hua X, Sun S, Mao Y, Song Y, Wang J, Li S. Induced CD10 expression during monocyte-to-macrophage differentiation identifies a unique subset of macrophages in pancreatic ductal adenocarcinoma. *Biochem Biophys Res Commun*. 2020; 524(4): 1064-1071. <https://doi.org/10.1016/j.bbrc.2020.02.042>.
11. Ino K, Suzuki T, Uehara C, Nagasaka T, Okamoto T, Kikkawa F, Mizutani S. The expression and localization of neutral endopeptidase 24.11/CD10 in human gestational trophoblastic diseases. *Lab Invest*. 2000; 80(11): 1729-38. <https://doi.org/10.1038/labinvest.3780183>.
12. Cozzolino I, Varone V, Picardi M, Baldi C, Memoli D, Ciancia G, Selleri C, De Rosa G, Vetrani A, Zeppa P. CD10, BCL6, and MUM1 expression in diffuse large B-cell lymphoma on FNA samples. *Cancer Cytopathol*. 2016; 124(2): 135-43. <https://doi.org/10.1002/cncy.21626>.
13. Hans CP, Weisenburger DD, Greiner TC, Gascoyne RD, Delabie J, Ott G, Muller-Hermelink HK, Campo E, Braziel RM,



- Jaffe ES, Pan Z, Farinha P, Smith LM, Falini B, Banham AH, Rosenwald A, Staudt LM, Connors JM, Armitage JO, Chan WC. Confirmation of the molecular classification of diffuse large B-cell lymphoma by immunohistochemistry using a tissue microarray. *Blood*. 2004; 103(1): 275-82. <https://doi.org/10.1182/blood-2003-05-1545>.
14. Tyagi A, Abrari A, Khurana A, Tyagi S. Immunohistochemical subtyping of diffuse large B-cell lymphoma into germinal center B-cell and activated B-cell subtype, along with correlation of the subtypes with extranodal involvement, serum lactate dehydrogenase, and positron emission tomography scan-based response assessment to chemotherapy. *J Cancer Res Ther*. 2022; 18(4): 1129-1136. https://doi.org/10.4103/jcrt.JCRT_842_20.
 15. Almasri NM, Iturraspe JA, Braylan RC. CD10 expression in follicular lymphoma and large cell lymphoma is different from that of reactive lymph node follicles. *Arch Pathol Lab Med*. 1998; 122(6): 539-44.
 16. Xu Y, McKenna RW, Kroft SH. Assessment of CD10 in the diagnosis of small B-cell lymphomas: a multiparameter flow cytometric study. *Am J Clin Pathol*. 2002; 117(2): 291-300. <https://doi.org/10.1309/T88X-71U4-WC0R-2531>.
 17. Chuang SS, Ye H, Du MQ, Lu CL, Dogan A, Hsieh PP, Huang WT, Jung YC. Histopathology and immunohistochemistry in distinguishing Burkitt lymphoma from diffuse large B-cell lymphoma with very high proliferation index and with or without a starry-sky pattern: a comparative study with EBER and FISH. *Am J Clin Pathol*. 2007; 128(4): 558-64. <https://doi.org/10.1309/EQJR3D3V0CCQGP04>.
 18. Yasir S, Herrera L, Gomez-Fernandez C, Reis IM, Umar S, Leveille R, Kava B, Jorda M. CD10+ and CK7/ROn-immunophenotype distinguishes renal cell carcinoma, conventional type with eosinophilic morphology from its mimickers. *Appl Immunohistochem Mol Morphol*. 2012; 20(5): 454-61. <https://doi.org/10.1097/PAI.0b013e31823fecd3>.
 19. Chu PG, Ishizawa S, Wu E, Weiss LM. Hepatocyte antigen as a marker of hepatocellular carcinoma: an immunohistochemical comparison to carcinoembryonic antigen, CD10, and alpha-fetoprotein. *Am J Surg Pathol*. 2002; 26(8): 978-88. <https://doi.org/10.1097/00000478-200208000-00002>.
 20. Agarwal A, Handa U, Kundu R, Sachdev A, Kochhar S. Hepatocyte paraffin-1, CD10, and CD34 immunostaining as a diagnostic aid in cytologic diagnosis of hepatic cancer. *J Cancer Res Ther*. 2022; 18(Supplement): S434-S438. https://doi.org/10.4103/jcrt.JCRT_467_20.
 21. Fleischmann A, Schlomm T, Huland H, Kollermann J, Simon P, Mirlacher M, Salomon G, Chun FH, Steuber T, Simon R, Sauter G, Graefen M, Erbersdobler A. Distinct subcellular expression patterns of neutral endopeptidase (CD10) in prostate cancer predict diverging clinical courses in surgically treated patients. *Clin Cancer Res*. 2008; 14(23): 7838-42. <https://doi.org/10.1158/1078-0432.CCR-08-1432>.
 22. Rapini RP, Special Stains, in *Practical dermatopathology*. 2021, Elsevier Health Sciences. p. 433-444.
 23. Lott R, Tunnicliffe J, Sheppard E, Santiago J, Hladik C, Nasim M, Zeitner K, Haas T, Kohl S, Movahedi-Lankarani S. *Practical Guide to Specimen Handling in Surgical Pathology*. 2022. College of American Pathologists (CAP) and National Society for Histotechnology (NSH): Northfield, IL.











Contact



If located within the United States, contact Sakura Finetek USA, Inc. by calling toll free **1-800-725-8723** or contact your Sakura Finetek representative or authorized distributor.

In countries other than the United States, contact the nearest authorized Sakura Finetek instrument distributor or representative. Contact details may be found at www.sakura.com

Any incident should be reported to the manufacturer. In the European Union, any serious incident that has occurred in relation to the device shall be reported to the manufacturer, authorized representative, and the competent authority of the appropriate Member State in which the user and/or the patient is established.

Symbols

-  Catalog number
-  Batch code
-  *in vitro* diagnostic medical device
-  Temperature limitation
-  Use by
-  Manufacturer
-  Consult instructions for use
-  European Conformity
-  Authorized representative in the European Community
-  Do not re-use (applies to capsule format only)

  Please see product label for lot and expiration date information and if available the date of manufacture

Storage: 2°C  8°C



CE

0123

	Sakura Finetek USA, Inc. 1750 W 214 th Street Torrance, CA 90501 U.S.A.
	Sakura Finetek Europe B.V. Flemingweg 10a 2408 AV Alphen aan den Rijn The Netherlands
Made in U.S.A.	

GS-33590 Rev. B, 2024-11
TE-13-031-01 v1.2

