

#### IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING **SECTION 1**

1.1. Product identifier

: TISSUE-TEK GENIE PRO DETECTION KIT, DAB DAB INTENSIFIER Product name

Product code : 8826-K250 intensifie

1.2. Relevant identified uses of the substance or mixture and uses advised against

: SU22 Professional use. For industrial or institutional use. PC21 Laboratory chemicals. Application

1.3. Details of the supplier of the safety data sheet

Supplier : Sakura Finetek Europe B.V.

Flemingweg 10A

2408 AV Alphen aan den Rijn, The Netherlands

Telephone +31-88-5920000 E-mail Application@sakura.eu Website www.sakura.eu

1.4. Emergency telephone number

EMERGENCY TELEPHONE NUMBER, for DOCTORS/FIRE BRIGADE/POLICE only:

NL - Telephone : +31-88-5920000 (During office hours only)

#### **SECTION 2** HAZARDS IDENTIFICATION

# 2.1. Classification of the substance or mixture

CLP classification : Not classified as dangerous according to Regulation (EC) No 1272/2008.

(1272/2008/EC)

Human health hazards : Low hazard for usual industrial or commercial handling.

Physical/chemical hazards : Not classified as dangerous according to statutory EC-Directives. Environmental hazards : Not classified as dangerous according to statutory EC-Directives.

2.2. Label elements

Label elements (1272/2008/EC):

Hazard pictograms : None.

Signal word : Not applicable.

H- and P-phrases : None.

2.3. Other hazards

Other information : Does not contain PBT or vPvB substances in concentrations higher than 0,1%.

#### **SECTION 3 COMPOSITION / INFORMATION ON INGREDIENTS**

3.2. Mixtures

Product description : Mixture.

Information on constituents : Does not contain hazardous substances in concentrations higher than 1%.

**SECTION 4** FIRST-AID MEASURES

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 1/7

Date of issue : 2019-01-25 Replaces issue dated INFO CARE SDS



# 4.1. Description of first aid measures

First aid measures

Inhalation : Move victim into fresh air. Consult a doctor if victim feels unwell.

Skin contact : Take off contaminated clothing. Wash off skin with plenty of water and soap before product dries up. Eye contact : Wash out with (lukewarm) water. Remove contact lenses. Consult a doctor if irritation persists.

Ingestion : Do not induce vomiting. Do rinse the mouth. Give one glass of water. Never give anything by mouth to an

unconscious person. Consult a doctor if victim feels unwell.

## 4.2. Most important symptoms and effects, both acute and delayed

Effects and symptoms

Inhalation : No specific effects and/or symptoms are known. Skin contact : No specific effects and/or symptoms are known.

Eye contact : May cause stinging of eyes and redness.

Ingestion : May cause a feeling of sickness, vomiting and diarrhoea.

## 4.3. Indication of any immediate medical attention and special treatment needed

Note to physicians : None known.

# SECTION 5 FIRE-FIGHTING MEASURES

## 5.1. Extinguishing media

Extinguishing media

Suitable : Carbondioxide (CO2). Foam. Dry chemical. Water fog.

Not suitable : None known.

# 5.2. Special hazards arising from the substance or mixture

Special exposure hazards : None known.

Hazardous thermal : Carbon monoxide may be evolved if incomplete combustion occurs.

decomposition products

## 5.3. Advice for firefighters

Special protective equipment : Use adequate respiratory equipment in case of insufficient ventilation.

for fire-fighters

# SECTION 6 ACCIDENTAL RELEASE MEASURES

# 6.1. Personal precautions, protective equipment and emergency procedures

Personal precautions : Danger of slipping. Clean up spills immediately. Wear shoes with non-slip soles.

# 6.2. Environmental precautions

Environmental precautions : Avoid release of product into sewers, surface water and/or ground water. In case of large spills: contain with

dike.

# 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Collect spilled material in containers. Absorb residues in sand or other inert material. Dispose at an authorised

waste collection point. Wash away remainder with plenty of water and soap.

# 6.4. Reference to other sections

Reference to other sections : See also section 8.

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 2/7

Date of issue : 2019-01-25 Replaces issue dated : --- INFO CARE SDS



#### **SECTION 7** HANDLING AND STORAGE

## 7.1. Precautions for safe handling

: Handle in accordance with good occupational hygiene and safety practices in well-ventilated areas. Avoid Handling

contact with skin and eyes.

## 7.2. Conditions for safe storage, including any incompatibilities

: Keep frost-free, in a cool, dry and well-ventilated place (< 35 °C).

Recommended packaging : Keep only in the original container. Non recommended packaging : Steel (except stainless steel).

7.3. Specific end use(s)

Use : Use only as directed.

#### **SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION**

## 8.1. Control parameters

Occupational exposure limits : Occupational exposure limits have not been established for this product. Derived no-effect levels (DNEL)

have not been established for this product. Predicted no-effect concentrations (PNEC) have not been

established for this product.

## 8.2. Exposure controls

Engineering measures : Use only in well-ventilated areas. Comply with standard precautionary measures for working with chemicals.

Hygienic measures When using do not eat, drink or smoke.

Personal protective equipment:

The efficiency of personal protective equipment depends among other things on temperature and degree of ventilation. Always get professional advice for the particular local situation.

Body protection : Use of specific protective industrial clothing is not required under normal conditions of use.

Respiratory protection Take care of sufficient ventilation.

Hand protection Under normal conditions of use specific gloves are not required. Wear appropriate gloves in case of frequent

or prolonged use and in case of large scale exposure. Suitable material: laminated film. ± 0.5 mm. Indication

of permeation breakthrough time: 6 hours.

Eye protection : Wear appropriate safety glasses when there is danger of possible eye contact.

#### **SECTION 9** PHYSICAL AND CHEMICAL PROPERTIES

## 9.1. Information on basic physical and chemical properties

Appearance : Liquid. Colour Blue

Odour Characteristic. Odour threshold Not known. : 2 - 11,5 рН : Soluble. Solubility in water

Partition coefficient (n-oc-

· Not known

Not measured. Not relevant for mixtures.

tanol/water)

Flash point : Not applicable. Does not contain combustible substances.

Flammability (solid, gas) : Not applicable. Liquid.

: Not known. Auto ignition temperature Does not contain substances with a known auto ignition temperature.

: > 90 °C Boiling point/boiling range Melting point/melting range : Not known.

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 3/7

Date of issue : 2019-01-25 Replaces issue dated INFO CARE SDS



The solvent content of this product is less than 1%.

Explosive properties : None known. Does not contain explosives.

Explosion limits (% in air) : Not applicable.

Oxidising properties Not applicable. Does not contain oxidizing substances.

Decomposition temperature Not known. Viscosity (20°C) Not known.

Viscosity (40°C) : Not relevant. The product contains < 10% substances having an aspiration hazard.

Vapour pressure (20°C) : Not known. Vapour density (20°C) Not relevant. Relative density (20°C) : 1 g/ml

Evaporation rate : Not known. (n-butyl acetate = 1)

9.2. Other information

Other information : Not relevant.

#### **SECTION 10** STABILITY AND REACTIVITY

10.1. Reactivity

Reactivity : See sub-sections below.

10.2. Chemical stability

Stability : Stable under normal conditions.

10.3. Possibility of hazardous reactions

Reactivity : No other hazardous reactions known.

10.4. Conditions to avoid

Conditions to avoid : See section 7.

10.5. Incompatible materials

Materials to avoid : No specific recommendations.

10.6. Hazardous decomposition products

Hazardous decomposition : Not known.

products

#### TOXICOLOGICAL INFORMATION **SECTION 11**

# 11.1. Information on toxicological effects

No toxicological research has been carried out on this product.

Inhalation

Acute toxicity : Calculated LC50: no data mg/l. Ingredients of unknown toxicity: < 1 %. ATE: > 5 mg/l. Low toxicity. Not

classified - based on available data, the classification criteria are not met. No specific effects and/or symptoms

are known.

Corrosion/irritation : No specific effects and/or symptoms are known. Not classified - based on available data, the classification

criteria are not met.

Sensitisation Does not contain substances classified as respiratory sensitiser. Not classified - based on available data, the

classification criteria are not met.

Carcinogenicity Not classified - based on available data, the classification criteria are not met. Mutagenicity : Not classified - based on available data, the classification criteria are not met.

Skin contact

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 4/7

Date of issue : 2019-01-25 Replaces issue dated INFO CARE SDS



Acute toxicity : Calculated LD50: > 5000 mg/kg.bw. Ingredients of unknown toxicity: < 1 %. ATE: > 2000 mg/kg.bw. Low

toxicity. Not classified - based on available data, the classification criteria are not met.

Corrosion/irritation : No specific effects and/or symptoms are known. Not classified - based on available data, the classification

criteria are not met.

Sensitisation : Does not contain skin sensitisers. Not classified - based on available data, the classification criteria are not

met

Mutagenicity : Does not contain mutagenic substances. Not classified - based on available data, the classification criteria are

not met

Eye contact

Corrosion/irritation : Slight irritation possible. Not classified - based on available data, the classification criteria are not met.

Ingestion

Acute toxicity : Calculated LD50: > 5000 mg/kg.bw. Ingredients of unknown toxicity: < 1 %. ATE: > 2000 mg/kg.bw. Low

toxicity. Not classified - based on available data, the classification criteria are not met.

Aspiration : Not classified - based on available data, the classification criteria are not met. Does not contain substances

with an aspiration hazard.

Corrosion/irritation : May cause a feeling of sickness, vomiting and diarrhoea.

Carcinogenicity : Does not contain carcinogenic substances. Not classified - based on available data, the classification criteria

are not met.

Mutagenicity : Does not contain mutagenic substances. Not classified - based on available data, the classification criteria are

not met.

Reprotoxicity : Development: Not expected to be reprotoxic. Development: Not classified - Based on available data, the

classification criteria are not met. Fertility: not expected to be reprotoxic. Fertility: Not classified - based on

available data, the classification criteria are not met.

# SECTION 12 ECOLOGICAL INFORMATION

# 12.1. Toxicity

No ecotoxicological research has been carried out on this product.

Ecotoxicity : Calculated LC50 (fish): 6 mg/l. Calculated EC50 (waterflea): 5050 mg/l. Contains 0 % of components with

unknown hazards to the aquatic environment. Not classified - based on available data, the classification

criteria are not met.

# 12.2. Persistence and degradability

Persistence – degradability : No specific information known. The surfactants contained in this preparation do not comply with the

biodegradability criteria as laid down in Regulation (EC) 648/2004 on detergents. This product shall not be

used in cleaners.

## 12.3. Bioaccumulative potential

Bioaccumulative potential : No specific information known.

12.4. Mobility in soil

Mobility : If product enters soil, it will be highly mobile and may contaminate groundwater.

# 12.5. Results of PBT and vPvB ass

PBT/vPvB assessment : Does not contain PBT or vPvB substances in concentrations higher than 0,1%.

12.6. Other adverse effects

Other information : Not applicable.

# SECTION 13 DISPOSAL CONSIDERATIONS

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 5/7

Date of issue : 2019-01-25 Replaces issue dated : --- INFO CARE SDS



### 13.1. Waste treatment methods

Product residues : Do not dispose empty pack with waste produced by households. Containers may be recycled. Treat product

residues and non-empty pack as chemical waste. Dispose waste to an official chemical waste depot.

Additional warning : None

Waste water discharge : Do not dispose into the environment, in drains or in water courses.

Local legislation : Disposal should be in accordance with applicable regional, national, and local laws and regulations. Local

regulations may be more stringent than regional or national requirements and must be complied with.

## SECTION 14 TRANSPORT INFORMATION

14.1. UN number

UN nr. : None.

14.2. UN proper shipping name

Transport name : Not regulated.

14.3/14.4/14.5. Transport hazard class(es)/Packing group/Environmental hazards

ADR/RID/ADN (road/railway/inland waterways)

Class : This product is not classified according to ADR/RID/ADN.

IMDG (sea)

Class : This product is not classified according to IMDG.

Marine pollutant : No

IATA (air)

Class : This product is not classified according to IATA.

14.6. Special precautions for user

Other information : Country specific variations may apply.

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Marpol : Not intended to be carried in bulk according to International Maritime Organisation (IMO) instruments.

Packaged liquids are not considered bulk.

## SECTION 15 REGULATORY INFORMATION

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Community regulations : Regulation (EU) No 2015/830 (REACH), Regulation (EC) No 1272/2008 (CLP) and other regulations.

15.2. Chemical safety assessment

Chemical safety assessment : Not applicable.

# SECTION 16 OTHER INFORMATION

## 16.1. Other information

The information in this safety data sheet is compiled in compliance with Regulation (EU) No 2015/830 dated 28 May 2015 and accurate to the best of our knowledge and experience at the date of issue specified. It is the user's obligation to use this product safely and to comply with all applicable laws and regulations concerning the use of the product. This safety data sheet complements the technical information sheets but does not replace them and offers no warranty with regard to product properties.

Users are also forewarned for any hazards involved when the product is used for other purposes than those for which it is designed.

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 6/7

Date of issue : 2019-01-25 Replaces issue dated : --- INFO CARE SDS



Changed or new information with regard to the previous release is indicated with an asterisk (\*).

List of abbreviations and acronyms that could be (but not necessarily are) used in this safety data sheet:

ADR : European Agreement concerning the International Carriage of Dangerous Goods by Road

ATE : Acute Toxicity Estimate

CLP : Classification, Labeling & Packaging

CMR : Carcinogenic, Mutagenic or toxic for Reproduction

EEC : European Economic Community

GHS : Globally Harmonized System of Classification and Labelling of Chemicals

IATA : International Air Transport Association IBC code : International Bulk Chemical Code

IMDG : International Maritime Dangerous Goods Code LD50/LC50 : Lethal Dose/Concentration for 50% of a population

MAC : Maximum Allowable Concentration

MARPOL : International Convention for the Prevention of Pollution From Ships

NO(A)EL : No Observed (Adverse) Effect Level

OECD : Organisation for Economic Co-operation and Development

PBT : Persistent, Bioaccumulative and Toxic

PC : Chemical product category

PT : Product type

REACH : Registration, Evaluation, Authorisation and Restriction of Chemicals

RID : Regulations concerning the International Carriage of Dangerous Goods by Rail

STP : Sewage Treatment Plant

SU : Sector of Use

TWA/STEL : Time-Weighted Average/Short Term Exposure Limit

UN : United Nations

VOC : Volatile Organic Compounds

vPvB : Very Persistent and Very Bioaccumulative

Procedure used to derive the classification according to Regulation (EC) No. 1272/2008:

Not classified : Based on test methods, experts judgement, bridging principles and calculation methods.

Full text of hazard classes mentioned in section 3: Not applicable.

Full text of H-phrases mentioned in section 3: Not applicable.

Number format : "," used as decimal separator.

End of safety data sheet.

Product name : Tissue-Tek Genie Pro Detection Kit, DAB DAB Intensifier Page 7/7

Date of issue : 2019-01-25 Replaces issue dated : --- INFO CARE SDS